

Girbau MS 613/617/623 LC Washer-extractor

Programme	Pre-wash	Main wash	Wash action	Number of rinses	Spin
Pr 1 Healthcare	Not heated	65°C	vigorous	3 Rinses	7 mins
Pr 2 Heavy soil-white	30°C	60°C	vigorous	4 Rinses	7 mins
Pr 3 Medium soil-white	30°C	50°C	vigorous	4 Rinses	7 mins
Pr 4 Light soil-white	No prewash	40°C	vigorous	4 Rinses	7 mins
Pr 5 Heavy soil – Colour	30°C	40°C	Normal	3 Rinses	7 mins
Pr 6 Medium soil – Colour	30°C	40°C	Normal	3 Rinses	7 mins
Pr 7 Synthetics	No prewash	40°C	Normal	3 Rinses	6 mins
Pr 8 Wool & delicates	No prewash	30°C	Delicate	3 Rinses	3 mins